

Product Overview

- Adjustable height from 71 cm to 123 cm (28 to 48.4 inches)
- 2 liftable columns - LCD height display with 3x heights memories
- Worksurface from W 100 to 190 x D80 cm (W 39.4 to 74.8 x D 31.5 inches)
- Weight capacity: 90 kg
- Worksurface can easily be lifted up or down by commuting an electrical switch
- Recommended use in the mass spectrometry field & the microscopy field
- Caster wheels allow to move the IonDesk easily (available with standard feet also)
- AC 110-220 V – 50-60 Hz



Work easily

in standing position...



...or in sitting position !



Description	Part-number
IonDesk 3 - 100x80	ID3-100
IonDesk 3 - 110x80	ID3-110
IonDesk 3 - 120x80	ID3-120
IonDesk 3 - 130x80	ID3-130
IonDesk 3 - 140x80	ID3-140
IonDesk 3 - 150x80	ID3-150
IonDesk 3 - 160x80	ID3-160
IonDesk 3 - 170x80	ID3-170
IonDesk 3 - 180x80	ID3-180
IonDesk 3 - 190x80	ID3-190
Optional Accessories	Part-number
CPU Holder	PC1
Optional Information Control Package (Flat screen monitor arm, keyboard & mouse platform)	ARM1

To be assembled